

## Product datasheet for **TL509236V**

### Wdr45 Mouse shRNA Lentiviral Particle (Locus ID 54636)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Wdr45 Mouse shRNA Lentiviral Particle (Locus ID 54636)
Locus ID:	54636
Synonyms:	C79260; DXlmx38e; JM5; Sfc19; Wdrx1; WIPI-4
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Wdr45 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC011479</a> , <a href="#">NM_001290792</a> , <a href="#">NM_001290794</a> , <a href="#">NM_001290795</a> , <a href="#">NM_172372</a> , <a href="#">NM_172372.1</a> , <a href="#">NM_172372.2</a> , <a href="#">NM_001290795.1</a> , <a href="#">NM_001290792.1</a> , <a href="#">NM_001290794.1</a>
UniProt ID:	<a href="#">Q91VM3</a>
Summary:	Component of the autophagy machinery that controls the major intracellular degradation process by which cytoplasmic materials are packaged into autophagosomes and delivered to lysosomes for degradation. Activated by the STK11/AMPK signaling pathway upon starvation, WDR45 is involved in autophagosome assembly downstream of WIPI2, regulating the size of forming autophagosomes. Probably recruited to membranes through its PtdIns3P activity. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).