

Product datasheet for **TL509495**

Avp Mouse shRNA Plasmid (Locus ID 11998)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Avp Mouse shRNA Plasmid (Locus ID 11998)
Locus ID:	11998
Synonyms:	Vp; Vsp
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Avp - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11998). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC051997 , NM_009732 , NM_009732.1 , NM_009732.2
UniProt ID:	P35455
Summary:	This gene encodes a member of the vasopressin/oxytocin family and preproprotein that is proteolytically processed to generate multiple protein products. These products include the neuropeptide hormone arginine vasopressin, and two other peptides, neurophysin 2 and copeptin. Arginine vasopressin binds to vasopressin receptors and functions as a vasopressor, to constrict blood vessels and increase blood pressure, and as an antidiuretic, to reduce the production of urine. Neurophysin 2 functions as a carrier protein in the transport of arginine vasopressin. This gene is present in a gene cluster with the related gene oxytocin on chromosome 2. [provided by RefSeq, Aug 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).