

## Product datasheet for **TL509671**

### Dsg2 Mouse shRNA Plasmid (Locus ID 13511)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Dsg2 Mouse shRNA Plasmid (Locus ID 13511)
Locus ID:	13511
Synonyms:	AA408168; D18Ertd293; D18Ertd293e
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Dsg2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 13511). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC131653</a> , <a href="#">BC131654</a> , <a href="#">NM_007883</a> , <a href="#">NM_007883.1</a> , <a href="#">NM_007883.3</a> , <a href="#">BC020144</a> , <a href="#">BC034056</a> , <a href="#">BC044848</a>
UniProt ID:	<a href="#">O55111</a>
Summary:	This gene encodes a member of the cadherin family of proteins that forms an integral transmembrane component of desmosomes, the multiprotein complexes involved in cell adhesion, organization of the cytoskeleton, cell sorting and cell signaling. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional protein. Mice lacking the encoded protein die in utero. Mutant mice lacking a part of the extracellular adhesive domain of the encoded protein develop cardiac fibrosis and dilation. This gene is located in a cluster of desmosomal cadherin genes on chromosome 18. [provided by RefSeq, Jan 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).