

Product datasheet for **TL510190**

Adamts12 Mouse shRNA Plasmid (Locus ID 239337)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Adamts12 Mouse shRNA Plasmid (Locus ID 239337)
Locus ID:	239337
Synonyms:	ADAMTS-12; AI605170
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Adamts12 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 239337). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_175501 , NM_175501.1 , NM_175501.2 , NM_175501.3 , BC151248 , BC152338
UniProt ID:	Q811B3
Summary:	This gene encodes a member of "a disintegrin and metalloproteinase with thrombospondin motifs" (ADAMTS) family of multi-domain matrix-associated metalloendopeptidases that have diverse roles in tissue morphogenesis and pathophysiological remodeling, in inflammation and in vascular biology. The encoded preproprotein undergoes proteolytic processing to generate an active protease. Mice lacking the encoded protein exhibit increased angiogenic response and tumor invasion in different models of angiogenesis and, severe inflammation and delayed recovery when subjected to experimental conditions that induce colitis, endotoxic sepsis and pancreatitis. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).