

Product datasheet for **TL510202**

Dxo Mouse shRNA Plasmid (Locus ID 112403)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Dxo Mouse shRNA Plasmid (Locus ID 112403)
Locus ID:	112403
Synonyms:	Dom3z; NG6
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Dxo - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 112403). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC004713 , NM_001163770 , NM_033613 , NM_033613.1 , NM_033613.2 , NM_001163770.1
UniProt ID:	O70348
Summary:	Ribonuclease that specifically degrades pre-mRNAs with a defective 5' end cap and is part of a pre-mRNA capping quality control. Has decapping, pyrophosphohydrolase and 5'-3' exonuclease activities. Has decapping activity toward incomplete 5' end cap mRNAs such as unmethylated 5' end-capped RNA to release GpppN and 5' end monophosphate RNA. The 5' end monophosphate RNA is then degraded by the 5'-3' exoribonuclease activity, enabling this enzyme to decap and degrade incompletely capped mRNAs. Also possesses RNA 5'-pyrophosphohydrolase activity by hydrolyzing the 5' end triphosphate to release pyrophosphates.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).