

Product datasheet for **TL510295**

Primpol Mouse shRNA Plasmid (Locus ID 408022)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Primpol Mouse shRNA Plasmid (Locus ID 408022)
Locus ID:	408022
Synonyms:	BC065112; Ccdc111
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Primpol - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 408022). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC065112 , NM_001001184 , NR_151680 , NM_001001184.1 , NM_001001184.2
UniProt ID:	Q6P1E7
Summary:	DNA primase and DNA polymerase able to initiate de novo DNA synthesis using dNTPs. Shows a high capacity to tolerate DNA damage lesions such as 8oxoG and abasic sites in DNA. Involved in translesion synthesis via its primase activity by mediating uninterrupted fork progression after programmed or damage-induced fork arrest and by reinitiating DNA synthesis after dNTP depletion. Required for mitochondrial DNA (mtDNA) synthesis, suggesting it may be involved in DNA tolerance during the replication of mitochondrial DNA. Has non-overlapping function with POLH (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).