

Product datasheet for **TL510438V**

Mttp Mouse shRNA Lentiviral Particle (Locus ID 17777)

Product data:

| | |
|---------------|---|
| Product Type: | shRNA Lentiviral Particles |
| Product Name: | Mttp Mouse shRNA Lentiviral Particle (Locus ID 17777) |
| Locus ID: | 17777 |
| Synonyms: | 1810043K16Rik; MTP |
| Vector: | pGFP-C-shLenti (TR30023) |
| Format: | Lentiviral particles |
| Components: | Mttp - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml. |
| RefSeq: | BC012686 , NM_001163457 , NM_008642 , NM_001355051 , NM_001355052 , NM_008642.1 , NM_008642.2 , NM_001163457.1 , NM_001163457.2 |
| UniProt ID: | O08601 |
| Summary: | Catalyzes the transport of triglyceride, cholesteryl ester, and phospholipid between phospholipid surfaces. Required for the secretion of plasma lipoproteins that contain apolipoprotein B. May be involved in regulating cholesteryl ester biosynthesis in cells that produce lipoproteins. Loads phospholipid into the C1D1 antigen-binding groove. Isoform 2 is critical for the development of natural killer T (NKT) cells. May have a role in the biogenesis of lipid droplets.[UniProtKB/Swiss-Prot Function] |
| shRNA Design: | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service . |

[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).