

## Product datasheet for **TL510461**

### Myo1g Mouse shRNA Plasmid (Locus ID 246177)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Myo1g Mouse shRNA Plasmid (Locus ID 246177)
Locus ID:	246177
Synonyms:	E430002D17Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Myo1g - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 246177). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC113755</a> , <a href="#">BC113954</a> , <a href="#">NM_178440</a> , <a href="#">NM_178440.1</a> , <a href="#">NM_178440.2</a> , <a href="#">NM_178440.3</a> , <a href="#">NM_178440.4</a> , <a href="#">BC028661</a>
UniProt ID:	<a href="#">Q5SUA5</a>
Summary:	Unconventional myosin required during immune response for detection of rare antigen-presenting cells by regulating T-cell migration (PubMed:25083865). Unconventional myosins are actin-based motor molecules with ATPase activity and serve in intracellular movements. Acts as a regulator of T-cell migration by generating membrane tension, enforcing cell-intrinsic meandering search, thereby enhancing detection of rare antigens during lymph-node surveillance, enabling pathogen eradication (PubMed:25083865). Also required in B-cells, where it regulates different membrane/cytoskeleton-dependent processes (PubMed:24310084). Involved in Fc-gamma receptor (Fc-gamma-R) phagocytosis (PubMed:23038771).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).