

## Product datasheet for **TL510714**

### Cacna1a Mouse shRNA Plasmid (Locus ID 12286)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cacna1a Mouse shRNA Plasmid (Locus ID 12286)
Locus ID:	12286
Synonyms:	alpha1A; APCA; BI; Caca1a; Cacn1a4; Cav2.1; Ccha1a; EA2; FHM; HPCA; Ia; MHP; MHP1; nmf352
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cacna1a - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 12286). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001252059</a> , <a href="#">NM_001252060</a> , <a href="#">NM_001252061</a> , <a href="#">NM_007578</a> , <a href="#">NM_007578.1</a> , <a href="#">NM_007578.2</a> , <a href="#">NM_007578.3</a> , <a href="#">NM_001252061.1</a> , <a href="#">NM_001252059.1</a> , <a href="#">NM_001252060.1</a>
UniProt ID:	<a href="#">P97445</a>
Summary:	Voltage-sensitive calcium channels (VSCC) mediate the entry of calcium ions into excitable cells and are also involved in a variety of calcium-dependent processes, including muscle contraction, hormone or neurotransmitter release, gene expression, cell motility, cell division and cell death. The isoform alpha-1A gives rise to P and/or Q-type calcium currents. P/Q-type calcium channels belong to the 'high-voltage activated' (HVA) group and are specifically blocked by the spider omega-agatoxin-IVA (AC P54282) (By similarity). They are however insensitive to dihydropyridines (DHP).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).