

## Product datasheet for **TL510717**

### Itgav Mouse shRNA Plasmid (Locus ID 16410)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Itgav Mouse shRNA Plasmid (Locus ID 16410)
Locus ID:	16410
Synonyms:	1110004F14Rik; 2610028E01Rik; CD51; D430040G12Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Itgav - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 16410). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_008402</a> , <a href="#">NM_008402.1</a> , <a href="#">NM_008402.2</a> , <a href="#">NM_008402.3</a> , <a href="#">BC021600</a> , <a href="#">BC048857</a> , <a href="#">BC133682</a> , <a href="#">BC167182</a>
UniProt ID:	<a href="#">P43406</a>
Summary:	This gene encodes a protein that is a member of the integrin superfamily. Integrins are transmembrane receptors involved cell adhesion and signaling, and they are subdivided based on the heterodimer formation of alpha and beta chains. This protein has been shown to heterodimerize with beta 1, beta 3, beta 6 and beta 8. The heterodimer of alpha v and beta 3 forms the Vitronectin receptor. This protein interacts with several extracellular matrix proteins to mediate cell adhesion and may play a role in cell migration. In mouse, deficiency of this gene is associated with defects in vascular morphogenesis in the brain and early post-natal death. [provided by RefSeq, May 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).