

## Product datasheet for **TL510719**

### Cdc42 Mouse shRNA Plasmid (Locus ID 12540)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cdc42 Mouse shRNA Plasmid (Locus ID 12540)
Locus ID:	12540
Synonyms:	AI747189; AU018915
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cdc42 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 12540). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC064792</a> , <a href="#">NM_001243769</a> , <a href="#">NM_009861</a> , <a href="#">NM_009861.1</a> , <a href="#">NM_009861.2</a> , <a href="#">NM_009861.3</a> , <a href="#">NM_001243769.1</a>
UniProt ID:	<a href="#">P60766</a>
Summary:	Plasma membrane-associated small GTPase which cycles between an active GTP-bound and an inactive GDP-bound state (PubMed:24352656). In its active state binds to a variety of effector proteins to regulate cellular responses. Involved in epithelial cell polarization processes. Regulates the bipolar attachment of spindle microtubules to kinetochores before chromosome congression in metaphase. Regulates cell migration (By similarity). In neurons, plays a role in the extension and maintenance of the formation of filopodia, thin and actin-rich surface projections. Required for DOCK10-mediated spine formation in Purkinje cells and hippocampal neurons (PubMed:25851601). Facilitates filopodia formation upon DOCK11-activation (PubMed:22494997). Upon activation by CaMKII, modulates dendritic spine structural plasticity by relaying CaMKII transient activation to synapse-specific, long-term signaling (By similarity). Also plays a role in phagocytosis through organization of the F-actin cytoskeleton associated with forming phagocytic cups (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).