

## Product datasheet for **TL510729**

### Pank4 Mouse shRNA Plasmid (Locus ID 269614)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Pank4 Mouse shRNA Plasmid (Locus ID 269614)
Locus ID:	269614
Synonyms:	D030031112Rik; R75150
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Pank4 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 269614). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC050089</a> , <a href="#">NM_001305804</a> , <a href="#">NM_172990</a> , <a href="#">NR_131230</a> , <a href="#">NM_172990.1</a> , <a href="#">NM_172990.2</a> , <a href="#">NM_172990.3</a> , <a href="#">NM_172990.4</a> , <a href="#">NM_172990.5</a> , <a href="#">BC050089.1</a> , <a href="#">BC012520</a> , <a href="#">BC019659</a> , <a href="#">BC027074</a>
UniProt ID:	<a href="#">Q80YV4</a>
Summary:	Plays a role in the physiological regulation of the intracellular CoA concentration (By similarity). The phosphatase activity shows preference for normal or oxidatively damaged intermediates of 4'-phosphopantetheine, which provides strong indirect evidence that the phosphatase activity pre-empts damage in the CoA pathway (By similarity). Hydrolyzing excess 4'-phosphopantetheine could constitute a directed overflow mechanism to prevent its oxidation to the S-sulfonate, sulfonate, or other forms (By similarity). Hydrolyzing 4'-phosphopantetheine sulfonate or S-sulfonate would forestall their conversion to inactive forms of CoA and acyl carrier protein (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).