

Product datasheet for **TL510922**

March5 Mouse shRNA Plasmid (Locus ID 69104)

Product data:

Product Type:	shRNA Plasmids
Product Name:	March5 Mouse shRNA Plasmid (Locus ID 69104)
Locus ID:	69104
Synonyms:	1810015H18Rik; 2310008I22Rik; 2700055A20Rik; 5730499H23Rik; E130202O05Rik; MARCH-V
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	March5 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 69104). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC064752 , BC107215 , BC132455 , BC132457 , NM_001164336 , NM_001164337 , NM_027314 , NM_027314.1 , NM_027314.2 , NM_027314.3 , NM_001164336.1 , NM_001164337.1 , BC107214
UniProt ID:	Q3KNM2
Summary:	Mitochondrial E3 ubiquitin-protein ligase that plays a crucial role in the control of mitochondrial morphology by acting as a positive regulator of mitochondrial fission. May play a role in the prevention of cell senescence acting as a regulator of mitochondrial quality control. Promotes ubiquitination of FIS1, DNM1L and MFN1.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).