

Product datasheet for **TL510958**

Mcu Mouse shRNA Plasmid (Locus ID 215999)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Mcu Mouse shRNA Plasmid (Locus ID 215999)
Locus ID:	215999
Synonyms:	2010012O16Rik; AV064928; C10orf42; Ccdc109a; D130073L02Rik; Gm64
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mcu - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 215999). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001033259 , NM_001033259.1 , NM_001033259.2 , NM_001033259.3 , BC139254 , BC139257
UniProt ID:	Q3UMR5



[View online »](#)

Summary: Mitochondrial inner membrane calcium uniporter that mediates calcium uptake into mitochondria (PubMed:21685886, PubMed:23900286, PubMed:24212091). Constitutes the pore-forming and calcium-conducting subunit of the uniporter complex (uni-plex) (By similarity). Activity is regulated by MICU1 and MICU2 (By similarity). At low Ca(2+) levels MCU activity is down-regulated by MICU1 and MICU2; at higher Ca(2+) levels MICU1 increases MCU activity (By similarity). Mitochondrial calcium homeostasis plays key roles in cellular physiology and regulates cell bioenergetics, cytoplasmic calcium signals and activation of cell death pathways (By similarity). Involved in buffering the amplitude of systolic calcium rises in cardiomyocytes (By similarity). While dispensable for baseline homeostatic cardiac function, acts as a key regulator of short-term mitochondrial calcium loading underlying a 'fight-or-flight' response during acute stress: acts by mediating a rapid increase of mitochondrial calcium in pacemaker cells (PubMed:26119742, PubMed:26119731, PubMed:25603276). Participates in mitochondrial permeability transition during ischemia-reperfusion injury (PubMed:26119731). Regulates glucose-dependent insulin secretion in pancreatic beta-cells by regulating mitochondrial calcium uptake (By similarity). Mitochondrial calcium uptake in skeletal muscle cells is involved in muscle size in adults (PubMed:25732818). Regulates synaptic vesicle endocytosis kinetics in central nerve terminal (PubMed:26644474). Involved in antigen processing and presentation (PubMed:25251370).[UniProtKB/Swiss-Prot Function]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).