

## Product datasheet for **TL511054**

### **Sgpl1 Mouse shRNA Plasmid (Locus ID 20397)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	Sgpl1 Mouse shRNA Plasmid (Locus ID 20397)
<b>Locus ID:</b>	20397
<b>Synonyms:</b>	AI428538; D10Xrf456; S1PL; Spl
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	Sgpl1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 20397). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">BC026135</a> , <a href="#">NM_001316673</a> , <a href="#">NM_001316674</a> , <a href="#">NM_009163</a> , <a href="#">NM_009163.1</a> , <a href="#">NM_009163.2</a> , <a href="#">NM_009163.3</a> , <a href="#">NM_009163.4</a>
<b>UniProt ID:</b>	<a href="#">Q8R0X7</a>
<b>Summary:</b>	Cleaves phosphorylated sphingoid bases (PSBs), such as sphingosine-1-phosphate, into fatty aldehydes and phosphoethanolamine (PubMed:9464245, PubMed:20097939). Elevates stress-induced ceramide production and apoptosis (PubMed:9464245). Required for global lipid homeostasis in liver and cholesterol homeostasis in fibroblasts (PubMed:20097939, PubMed:28262793). Involved in the regulation of pro-inflammatory response and neutrophil trafficking (PubMed:21173151). Modulates neuronal autophagy via phosphoethanolamine production which regulates accumulation of aggregate-prone proteins such as APP (PubMed:28521611). Seems to play a role in establishing neuronal contact sites and axonal maintenance (By similarity).[UniProtKB/Swiss-Prot Function]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).