

Product datasheet for **TL511168V**

Ube2i Mouse shRNA Lentiviral Particle (Locus ID 22196)

Product data:

| | |
|---------------|--|
| Product Type: | shRNA Lentiviral Particles |
| Product Name: | Ube2i Mouse shRNA Lentiviral Particle (Locus ID 22196) |
| Locus ID: | 22196 |
| Synonyms: | 5830467E05Rik; F830028O17Rik; UBC9; Ubce2i; Ubce9 |
| Vector: | pGFP-C-shLenti (TR30023) |
| Format: | Lentiviral particles |
| Components: | Ube2i - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml. |
| RefSeq: | BC037635 , NM_001177609 , NM_001177610 , NM_011665 , NM_011665.1 , NM_011665.2 , NM_011665.3 , NM_011665.4 , NM_001177610.1 , NM_001177609.1 , BC020437 , NM_001357694 |
| UniProt ID: | P63280 |
| Summary: | Accepts the ubiquitin-like proteins SUMO1, SUMO2 and SUMO3 from the UBLE1A-UBLE1B E1 complex and catalyzes their covalent attachment to other proteins with the help of an E3 ligase such as RANBP2, CBX4 and ZNF451. Can catalyze the formation of poly-SUMO chains. Essential for nuclear architecture, chromosome segregation and embryonic viability. Necessary for sumoylation of FOXL2 and KAT5 (By similarity). Sumoylates p53/TP53 at 'Lys-386'. [UniProtKB/Swiss-Prot Function] |
| shRNA Design: | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service . |




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL511168A virus into HEK293 cells. TL511168A virus was prepared using lenti-shRNA TL511168A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL511168B virus into HEK293 cells. TL511168B virus was prepared using lenti-shRNA TL511168B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL511168C] virus into HEK293 cells. [TL511168C] virus was prepared using lenti-shRNA [TL511168C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL511168D] virus into HEK293 cells. [TL511168D] virus was prepared using lenti-shRNA [TL511168D] and [TR30037] packaging kit.