

Product datasheet for **TL511180**

Fnip1 Mouse shRNA Plasmid (Locus ID 216742)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Fnip1 Mouse shRNA Plasmid (Locus ID 216742)
Locus ID:	216742
Synonyms:	A730024A03Rik; AI838773; AW557298
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Fnip1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 216742). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC079892 , NM_173753 , NM_173753.1 , NM_173753.2 , NM_173753.3 , NM_173753.4 , BC024776
UniProt ID:	Q68FD7
Summary:	Acts as a co-chaperone of HSP90AA1. Inhibits the ATPase activity of HSP90AA1 leading to reduction in its chaperone activity. Facilitates the binding of client protein FLCN to HSP90AA1. Competes with the activating co-chaperone AHSA1 for binding to HSP90AA1, thereby providing a reciprocal regulatory mechanism for chaperoning of client proteins. May be involved in energy and/or nutrient sensing through the AMPK and mTOR signaling pathways. May regulate phosphorylation of RPS6KB1.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).