

Product datasheet for **TL511239**

Prdm16 Mouse shRNA Plasmid (Locus ID 70673)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Prdm16 Mouse shRNA Plasmid (Locus ID 70673)
Locus ID:	70673
Synonyms:	5730557K01Rik; csp1; mel1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Prdm16 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 70673). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC059838 , NM_001177995 , NM_001291026 , NM_001291029 , NM_027504 , NM_001177995.1 , NM_027504.1 , NM_027504.2 , NM_027504.3 , NM_001291029.1 , NM_001291026.1 , BC053082
UniProt ID:	A2A935
Summary:	Binds DNA and functions as a transcriptional regulator (PubMed:18483224). Displays histone methyltransferase activity and monomethylates 'Lys-9' of histone H3 (H3K9me1) in vitro (PubMed:22939622). Probably catalyzes the monomethylation of free histone H3 in the cytoplasm which is then transported to the nucleus and incorporated into nucleosomes where SUV39H methyltransferases use it as a substrate to catalyze histone H3 'Lys-9' trimethylation (PubMed:22939622). Likely to be one of the primary histone methyltransferases along with MECOM/PRDM3 that direct cytoplasmic H3K9me1 methylation (PubMed:22939622). Functions in the differentiation of brown adipose tissue (BAT) which is specialized in dissipating chemical energy in the form of heat in response to cold or excess feeding while white adipose tissue (WAT) is specialized in the storage of excess energy and the control of systemic metabolism (PubMed:17618855, PubMed:18483224). Together with CEBPB, regulates the differentiation of myoblastic precursors into brown adipose cells (PubMed:18719582, PubMed:19641492). Functions as a repressor of TGF-beta signaling. [UniProtKB/Swiss-Prot Function]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).