

Product datasheet for **TL511309**

Krt16 Mouse shRNA Plasmid (Locus ID 16666)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Krt16 Mouse shRNA Plasmid (Locus ID 16666)
Locus ID:	16666
Synonyms:	AI324768; CK-16; K16; Krt1-1; Krt1-16
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Krt16 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 16666). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC103615 , BC103616 , BC103617 , BC103666 , BC114977 , NM_008470 , NM_008470.1
UniProt ID:	Q9Z2K1
Summary:	The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The encoded protein is a cytokeratin and acts as an innate immune system effector, promoting the inflammatory response upon breach of the skin barrier. Defects in this gene are a cause of pachyonychia congenita. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).