

Product datasheet for **TL511319**

Jmjd6 Mouse shRNA Plasmid (Locus ID 107817)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Jmjd6 Mouse shRNA Plasmid (Locus ID 107817)
Locus ID:	107817
Synonyms:	5730436I23Rik; D11Erttd195e; mKIAA0585; PSR; PtdSerR; Ptdsr
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Jmjd6 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 107817). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC056629 , NM_033398 , NM_033398.1 , NM_033398.2 , BC006067 , BC011503 , BC051491 , BC117737 , NM_001363363
UniProt ID:	Q9ER15
Summary:	This gene encodes a nuclear protein with a JmjC domain. JmjC domain-containing proteins are predicted to function as protein hydroxylases or histone demethylases. This protein functions in differentiation of multiple tissues during development, and in anti-inflammatory cytokine signaling. It was first identified as a putative phosphatidylserine receptor involved in phagocytosis of apoptotic cells; however, subsequent studies have indicated that this protein does not directly function in the clearance of apoptotic cells, and questioned whether it is a true phosphatidylserine receptor. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).