

Product datasheet for **TL511456**

Unc5c Mouse shRNA Plasmid (Locus ID 22253)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Unc5c Mouse shRNA Plasmid (Locus ID 22253)
Locus ID:	22253
Synonyms:	B130051O18Rik; rcm; Unc5h3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Unc5c - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 22253). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC115772 , NM_001293561 , NM_009472 , NM_009472.1 , NM_009472.2 , NM_009472.3 , NM_009472.4 , NM_001293561.1
UniProt ID:	O08747
Summary:	Receptor for netrin required for axon guidance (PubMed:22685302, PubMed:10399920). Mediates axon repulsion of neuronal growth cones in the developing nervous system upon ligand binding (PubMed:10399920, PubMed:22685302). NTN1/Netrin-1 binding might cause dissociation of UNC5C from polymerized TUBB3 in microtubules and thereby lead to increased microtubule dynamics and axon repulsion (PubMed:28483977). Axon repulsion in growth cones may also be caused by its association with DCC that may trigger signaling for repulsion (PubMed:10399920). Might also collaborate with DSCAM in NTN1-mediated axon repulsion independently of DCC (PubMed:22685302). Also involved in corticospinal tract axon guidance independently of DCC (PubMed:9126743, PubMed:9389662, PubMed:12451134). Involved in dorsal root ganglion axon projection towards the spinal cord (By similarity). It also acts as a dependence receptor required for apoptosis induction when not associated with netrin ligand (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).