

Product datasheet for **TL511780V**

Sox17 Mouse shRNA Lentiviral Particle (Locus ID 20671)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Sox17 Mouse shRNA Lentiviral Particle (Locus ID 20671)
Locus ID:	20671
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Sox17 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001289464 , NM_001289465 , NM_001289466 , NM_001289467 , NM_011441 , NM_011441.1 , NM_011441.2 , NM_011441.3 , NM_011441.4 , NM_011441.5 , NM_001289467.1 , NM_001289465.1 , NM_001289466.1 , NM_001289464.1 , BC027687 , BC037521 , BC060612 , BC148437 , BC153021
UniProt ID:	Q61473
Summary:	This gene encodes a member of the Sox (Sry-related high mobility group box) family of transcription factors involved in the regulation of embryonic development. The encoded protein plays a role in the determination of cell fate and in maintaining cell identity. This gene regulates tumor angiogenesis and tumor progression. Mutations in the human gene are associated with vesicoureteral reflux, characterized by the backward flow of urine from the bladder into the ureters or the kidney. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).