

Product datasheet for **TL512125**

Rab38 Mouse shRNA Plasmid (Locus ID 72433)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Rab38 Mouse shRNA Plasmid (Locus ID 72433)
Locus ID:	72433
Synonyms:	2310011F14Rik; AU043391; cht
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Rab38 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 72433). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC089578 , NM_028238 , NM_028238.1 , NM_028238.2 , NM_028238.3 , NM_028238.4 , NM_028238.5 , NM_028238.6 , NM_028238.7
UniProt ID:	Q8QZZ8
Summary:	Plays a role in the maturation of phagosomes that engulf pathogens, such as S.aureus and Mycobacterium (By similarity). May be involved in melanosomal transport and docking. Involved in the proper sorting of TYRP1. Involved in peripheral melanosomal distribution of TYRP1 in melanocytes; the function, which probably is implicating vesicle-trafficking, includes cooperation with ANKRD27 and VAMP7 (PubMed:21187289). Plays an important role in the control of melanin production and melanosome biogenesis (By similarity). In concert with RAB32, regulates the proper trafficking of melanogenic enzymes TYR, TYRP1 and DCT/TYRP2 to melanosomes in melanocytes (PubMed:26620560).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).