

## Product datasheet for **TL512189**

### Ints8 Mouse shRNA Plasmid (Locus ID 72656)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Ints8 Mouse shRNA Plasmid (Locus ID 72656)
Locus ID:	72656
Synonyms:	2810013E07Rik; AV063769; D130008D20Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Ints8 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 72656). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC049092</a> , <a href="#">BC050192</a> , <a href="#">BC057618</a> , <a href="#">NM_001159595</a> , <a href="#">NM_178112</a> , <a href="#">NM_001159595.1</a> , <a href="#">NM_178112.1</a> , <a href="#">NM_178112.2</a> , <a href="#">NM_178112.3</a> , <a href="#">NM_178112.4</a> , <a href="#">NM_178112.5</a> , <a href="#">BC024524</a> , <a href="#">BC060149</a>
UniProt ID:	<a href="#">Q80V86</a>
Summary:	Component of the Integrator complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3'-box-dependent processing. The Integrator complex is associated with the C-terminal domain (CTD) of RNA polymerase II largest subunit (POLR2A) and is recruited to the U1 and U2 snRNAs genes (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).