

Product datasheet for **TL512231**

Akt1 Mouse shRNA Plasmid (Locus ID 11651)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Akt1 Mouse shRNA Plasmid (Locus ID 11651)
Locus ID:	11651
Synonyms:	Ak; Akt; LTR-akt; PK; PKB; PKB/A; PKB/Akt; PKBalpha; Rac
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Akt1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11651). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC066018 , NM_001165894 , NM_009652 , NM_009652.1 , NM_009652.2 , NM_009652.3 , NM_001165894.1 , BC115583
UniProt ID:	P31750
Summary:	This gene encodes the founding member of the Akt serine-threonine protein kinase gene family that also includes Akt2 and Akt3. This kinase is a major downstream effector of the phosphatidylinositol 3-kinase (PI3K) pathway that mediates the effects of various growth factors such as platelet-derived growth factor (PDGF), epidermal growth factor (EGF), insulin and insulin-like growth factor I (IGF-I). It is activated through recruitment to cellular membranes by PI3K lipid products and by phosphorylation by 3-phosphoinositide dependent kinase-1. It then further phosphorylates different downstream proteins in response to various extracellular signals and thus plays a pivotal role in mediating a variety of cellular processes, such as glucose metabolism, glycogen biosynthesis, protein synthesis and turn over, inflammatory response, cell survival (anti-apoptosis) and development. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL512231A virus into HEK293 cells. TL512231A virus was prepared using lenti-shRNA TL512231A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL512231B virus into HEK293 cells. TL512231B virus was prepared using lenti-shRNA TL512231B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL512231C] virus into HEK293 cells. [TL512231C] virus was prepared using lenti-shRNA [TL512231C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL512231D] virus into HEK293 cells. [TL512231D] virus was prepared using lenti-shRNA [TL512231D] and [TR30037] packaging kit.