

## Product datasheet for **TL512258**

### Itga3 Mouse shRNA Plasmid (Locus ID 16400)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Itga3 Mouse shRNA Plasmid (Locus ID 16400)
Locus ID:	16400
Synonyms:	AA407068; CD49C; GAPB3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Itga3 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 16400). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC053031</a> , <a href="#">BC062205</a> , <a href="#">NM_013565</a> , <a href="#">NM_013565.1</a> , <a href="#">NM_013565.2</a> , <a href="#">NM_013565.3</a> , <a href="#">BC023259</a> , <a href="#">BC096741</a>
UniProt ID:	<a href="#">Q62470</a>
Summary:	This gene encodes a subunit of integrin family of cell surface proteins. The encoded protein undergoes post-translational processing to form a disulfide bond-linked dimer comprised of heavy and light chains. At the cell surface, the encoded protein non-covalently associates with the integrin beta-1 subunit to form a heterodimer that interacts with many extracellular matrix proteins including fibronectin and laminin. Mice lacking the encoded protein die during the first day after birth due to severe abnormalities in kidneys. Mice lacking the encoded protein specifically in the basal layer of epidermis display several skin defects and accelerated wound healing. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).