

Product datasheet for **TL512454V**

Cryaa Mouse shRNA Lentiviral Particle (Locus ID 12954)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cryaa Mouse shRNA Lentiviral Particle (Locus ID 12954)
Locus ID:	12954
Synonyms:	Acry; Acry-1; Cry; Crya; Crya-1; Crya1; DAcry; DAcry-1; lop18
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Cryaa - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC085170 , BC085172 , BC092385 , NM_001278569 , NM_001278570 , NM_013501 , NM_013501.1 , NM_013501.2 , NM_013501.3 , NM_001278570.1 , NM_001278569.1
UniProt ID:	P24622
Summary:	This gene encodes subunit a, one of two subunits of alpha-crystallin, which is a high molecular weight, soluble aggregate and is a member of the small heat shock protein (sHSP) family. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. It acts as a molecular chaperone and is the major protein in the eye lens, maintaining the transparency and refractive index of the lens. In mouse, deficiency in this gene is associated with smaller lenses and eyes and with increasing lens opacity with age. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).