

## Product datasheet for **TL512538**

### Appl1 Mouse shRNA Plasmid (Locus ID 72993)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Appl1 Mouse shRNA Plasmid (Locus ID 72993)
Locus ID:	72993
Synonyms:	2900057D21Rik; 7330406P05Rik; AI585782; AW209077; BB022931; C88264; DIP13
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Appl1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 72993). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC063751</a> , <a href="#">NM_145221</a> , <a href="#">NM_145221.1</a> , <a href="#">NM_145221.2</a> , <a href="#">BC019708</a>
UniProt ID:	<a href="#">Q8K3H0</a>
Summary:	Multifunctional adapter protein that binds to various membrane receptors, nuclear factors and signaling proteins to regulate many processes, such as cell proliferation, immune response, endosomal trafficking and cell metabolism (By similarity) (PubMed:25328665, PubMed:25568335, PubMed:27219021). Regulates signaling pathway leading to cell proliferation through interaction with RAB5A and subunits of the NuRD/MeCP1 complex (By similarity). Functions as a positive regulator of innate immune response via activation of AKT1 signaling pathway by forming a complex with APPL1 and PIK3R1 (PubMed:25328665). Inhibits Fc-gamma receptor-mediated phagocytosis through PI3K/Akt signaling in macrophages (PubMed:25568335). Regulates TLR4 signaling in activated macrophages (PubMed:27219021). Involved in trafficking of the TGFBR1 from the endosomes to the nucleus via microtubules in a TRAF6-dependent manner. Plays a role in cell metabolism by regulating adiponectin and insulin signaling pathways (By similarity). Required for fibroblast migration through HGF cell signaling (PubMed:26445298). Positive regulator of beta-catenin/TCF-dependent transcription through direct interaction with RUVBL2/reptin resulting in the relief of RUVBL2-mediated repression of beta-catenin/TCF target genes by modulating the interactions within the beta-catenin-reptin-HDAC complex (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).