

Product datasheet for **TL512547**

Rnf10 Mouse shRNA Plasmid (Locus ID 50849)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Rnf10 Mouse shRNA Plasmid (Locus ID 50849)
Locus ID:	50849
Synonyms:	AA675014; mKIAA0262; RIE2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Rnf10 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 50849). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC010342 , BC058527 , NM_001302448 , NM_001302449 , NM_016698 , NM_016698.1 , NM_016698.2 , NM_001302449.1 , NM_001302448.1
UniProt ID:	Q3UIW5
Summary:	The protein encoding this gene is a member of the really interesting new gene finger protein family. Members of this family contain protein motifs similar to zinc finger domains and are involved in many processes that include transcriptional regulation, DNA repair and signal transduction. Expression of this gene is upregulated during neuronal differentiation of cultured cells, and inhibition of its expression impairs differentiation and cell cycle exit, providing evidence for a function in neuronal differentiation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).