

## Product datasheet for **TL512561**

### Tmie Mouse shRNA Plasmid (Locus ID 20776)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Tmie Mouse shRNA Plasmid (Locus ID 20776)
Locus ID:	20776
Synonyms:	5131400L21Rik; Mm.87012; sr
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Tmie - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 20776). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_146260</a> , <a href="#">NM_146260.1</a> , <a href="#">NM_146260.2</a> , <a href="#">BC138840</a> , <a href="#">BC141124</a> , <a href="#">NM_146260.3</a>
UniProt ID:	<a href="#">Q8K467</a>
Summary:	Unknown. The protein may play some role in a cellular membrane location. May reside within an internal membrane compartment and function in pathways such as those involved in protein and/or vesicle trafficking. Alternatively, the mature protein may be localized in the plasma membrane and serve as a site of interaction for other molecules through its highly charged C-terminal domain.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).