

Product datasheet for **TL512579**

Tirap Mouse shRNA Plasmid (Locus ID 117149)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tirap Mouse shRNA Plasmid (Locus ID 117149)
Locus ID:	117149
Synonyms:	AA407980; C130027E04Rik; Mal; Tlr4ap; Wyatt
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Tirap - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 117149). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC005584 , NM_001177845 , NM_001177846 , NM_001177847 , NM_054096 , NM_054096.1 , NM_054096.2 , NM_001177845.1 , NM_001177846.1 , NM_001177847.1
UniProt ID:	Q99JY1
Summary:	Adapter involved in the TLR2 and TLR4 signaling pathways in the innate immune response. Acts via IRAK2 and TRAF-6, leading to the activation of NF-kappa-B, MAPK1, MAPK3 and JNK, and resulting in cytokine secretion and the inflammatory response (By similarity). Positively regulates the production of TNF-alpha and interleukin-6 (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).