

Product datasheet for **TL512758**

Tmem175 Mouse shRNA Plasmid (Locus ID 72392)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tmem175 Mouse shRNA Plasmid (Locus ID 72392)
Locus ID:	72392
Synonyms:	3010001K23Rik; AI504381; mTMEM175
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Tmem175 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 72392). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC005542 , NM_001163531 , NM_001163532 , NM_028223 , NM_001359423 , NM_001359426 , NM_001359427 , NM_001163532.1 , NM_028223.1 , NM_028223.2 , NM_028223.3 , NM_001163531.1
UniProt ID:	Q9CXY1
Summary:	Organelle-specific potassium channel specifically responsible for potassium conductance in endosomes and lysosomes. Forms a potassium-permeable leak-like channel, which regulates luminal pH stability and is required for autophagosome-lysosome fusion. Constitutes the major lysosomal potassium channel.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL512758A virus into HEK293 cells. TL512758A virus was prepared using lenti-shRNA TL512758A and [TR30037] packaging kit.




GFP signal was observed under microscope at 48 hours after transduction of TL512758A virus into HEK293 cells. TL512758A virus was prepared using lenti-shRNA TL512758A and [TR30037] packaging kit.




GFP signal was observed under microscope at 48 hours after transduction of TL512758B virus into HEK293 cells. TL512758B virus was prepared using lenti-shRNA TL512758B and [TR30037] packaging kit.




GFP signal was observed under microscope at 48 hours after transduction of TL512758B virus into HEK293 cells. TL512758B virus was prepared using lenti-shRNA TL512758B and [TR30037] packaging kit.




GFP signal was observed under microscope at 48 hours after transduction of [TL512758C] virus into HEK293 cells. [TL512758C] virus was prepared using lenti-shRNA [TL512758C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL512758C] virus into HEK293 cells. [TL512758C] virus was prepared using lenti-shRNA [TL512758C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL512758D] virus into HEK293 cells. [TL512758D] virus was prepared using lenti-shRNA [TL512758D] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL512758D] virus into HEK293 cells. [TL512758D] virus was prepared using lenti-shRNA [TL512758D] and [TR30037] packaging kit.