

## Product datasheet for **TL512784**

### **Klk13 Mouse shRNA Plasmid (Locus ID 626834)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	Klk13 Mouse shRNA Plasmid (Locus ID 626834)
<b>Locus ID:</b>	626834
<b>Synonyms:</b>	Egfbp-2; mGk-13
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	Klk13 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 626834). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">NM_001039042</a> , <a href="#">NM_001039042.1</a> , <a href="#">BC063763</a> , <a href="#">BC152778</a> , <a href="#">BC156599</a>
<b>Summary:</b>	Kallikreins are serine proteases encoded by a cluster of highly related genes on chromosome 7. When genomic sequence became available, it was possible to clarify the gene-to-sequence relationship for this family. Until January, 2006, NM_010115 was called kallikrein 13. After that time, NM_001039042 was called Klk13 and NM_010115 was called Klk1b26. [provided by RefSeq, Jul 2008]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).