

Product datasheet for **TL512813V**

Rbpj Mouse shRNA Lentiviral Particle (Locus ID 19664)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Rbpj Mouse shRNA Lentiviral Particle (Locus ID 19664)
Locus ID:	19664
Synonyms:	AI843960; CBF1; Igkjrj; Igkrsbp; RBP-J; RBP-J kappa; RBP-Jkappa; RBPjk; Rbpsuh
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Rbpj - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	<u>BC051387</u> , <u>NM_001080927</u> , <u>NM_001080928</u> , <u>NM_001277116</u> , <u>NM_009035</u> , <u>NM_001359152</u> , <u>NM_001080927.1</u> , <u>NM_001080927.2</u> , <u>NM_009035.1</u> , <u>NM_009035.3</u> , <u>NM_009035.4</u> , <u>NM_009035.5</u> , <u>NM_001080928.1</u> , <u>NM_001277116.1</u> , <u>BC035299</u>
UniProt ID:	<u>P31266</u>
Summary:	Transcriptional regulator that plays a central role in Notch signaling, a signaling pathway involved in cell-cell communication that regulates a broad spectrum of cell-fate determinations (PubMed:7566092). Acts as a transcriptional repressor when it is not associated with Notch proteins. When associated with some NICD product of Notch proteins (Notch intracellular domain), it acts as a transcriptional activator that activates transcription of Notch target genes. Probably represses or activates transcription via the recruitment of chromatin remodeling complexes containing histone deacetylase or histone acetylase proteins, respectively. Specifically binds to the immunoglobulin kappa-type J segment recombination signal sequence. Binds specifically to methylated DNA. Binds to the oxygen responsive element of COX4I2 and activates its transcription under hypoxia conditions (4% oxygen) (By similarity). Negatively regulates the phagocyte oxidative burst in response to bacterial infection by repressing transcription of NADPH oxidase subunits (PubMed:26194095).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).