

## Product datasheet for **TL512871**

### Mkks Mouse shRNA Plasmid (Locus ID 59030)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Mkks Mouse shRNA Plasmid (Locus ID 59030)
Locus ID:	59030
Synonyms:	1300013E18Rik; AI463362; AI957237; Bbs; Bbs6
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mkks - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 59030). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC117836</a> , <a href="#">NM_001141946</a> , <a href="#">NM_001286981</a> , <a href="#">NM_001286983</a> , <a href="#">NM_021527</a> , <a href="#">NM_001141946.1</a> , <a href="#">NM_021527.1</a> , <a href="#">NM_021527.2</a> , <a href="#">NM_001286983.1</a> , <a href="#">NM_001286981.1</a> , <a href="#">BC024359</a> , <a href="#">BC058174</a> , <a href="#">BC080765</a> , <a href="#">BC100336</a> , <a href="#">BC117835</a>
UniProt ID:	<a href="#">Q9J170</a>
Summary:	This gene encodes a protein which shares sequence similarity with other members of the type II chaperonin family. The encoded protein is a centrosome-shuttling protein and plays an important role in cytokinesis. This protein also interacts with other type II chaperonin members to form a complex known as the BBSome, which involves ciliary membrane biogenesis. This protein is encoded by a downstream open reading frame (dORF). Several upstream open reading frames (uORFs) have been identified, which repress the translation of the dORF, and two of which can encode small mitochondrial membrane proteins. Alternatively spliced transcripts encoding distinct isoforms have been found for this gene. [provided by RefSeq, Nov 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).