

## Product datasheet for **TL513058**

### Stard3 Mouse shRNA Plasmid (Locus ID 59045)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Stard3 Mouse shRNA Plasmid (Locus ID 59045)
Locus ID:	59045
Synonyms:	es64; Mln64
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Stard3 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 59045). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC003313</a> , <a href="#">NM_021547</a> , <a href="#">NM_021547.1</a> , <a href="#">NM_021547.2</a> , <a href="#">NM_021547.3</a> , <a href="#">BC058727</a> , <a href="#">BC114641</a>
UniProt ID:	<a href="#">Q61542</a>
Summary:	Sterol-binding protein that mediates cholesterol transport from the endoplasmic reticulum to endosomes (By similarity). Creates contact site between the endoplasmic reticulum and late endosomes: localizes to late endosome membranes and contacts the endoplasmic reticulum via interaction with VAPA and VAPB (By similarity). Acts as a lipid transfer protein that redirects sterol to the endosome at the expense of the cell membrane and favors membrane formation inside endosomes (By similarity). May also mediate cholesterol transport between other membranes, such as mitochondria membrane or cell membrane (By similarity). However, such results need additional experimental evidences; probably mainly mediates cholesterol transport from the endoplasmic reticulum to endosomes (By similarity). Does not activate transcriptional cholesterol sensing (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).