

Product datasheet for **TL513127**

Ano7 Mouse shRNA Plasmid (Locus ID 404545)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ano7 Mouse shRNA Plasmid (Locus ID 404545)
Locus ID:	404545
Synonyms:	IPC; NG; NGE; Ngep; Ngep-L; Pc; Pcanap5; Tmem16g
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Ano7 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 404545). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC116706 , NM_001271884 , NM_207031 , NM_207031.1 , NM_207031.2 , NM_001271884.1
UniProt ID:	Q14AT5
Summary:	This gene encodes a member of the anoctamin family, which in mammals is comprised of 10 members. Anoctamin proteins are proposed to have eight transmembrane domains with both termini facing the cytoplasm and a C-terminal domain of unknown function. While some members have been characterized as calcium-activated chloride channels, this protein is reported to have little anion conductance activity. In humans, this protein is primarily found in prostate tissues and may serve as a target for prostate cancer immunotherapy. Alternative splicing results in multiple transcript variants that encode different isoforms. [provided by RefSeq, Dec 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).