

Product datasheet for **TL513204**

Tfe3 Mouse shRNA Plasmid (Locus ID 209446)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tfe3 Mouse shRNA Plasmid (Locus ID 209446)
Locus ID:	209446
Synonyms:	bHLHe33; F830016E06Rik; Tcfe3; Tfe-3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Tfe3 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 209446). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC056358 , BC063047 , NM_001105196 , NM_001105197 , NM_001271489 , NM_001271490 , NM_001271491 , NM_172472 , NM_172472.1 , NM_172472.2 , NM_172472.3 , NM_001105196.1 , NM_001271491.1 , NM_001105197.1 , NM_001271489.1 , NM_001271490.1 , BC032963 , BC051568
UniProt ID:	Q64092
Summary:	Transcription factor that specifically recognizes and binds E-box sequences (5'-CANNTG-3') (PubMed:16936731). Efficient DNA-binding requires dimerization with itself or with another MiT/TFE family member such as TFEB or MITF (PubMed:16936731). In association with TFEB, activates the expression of CD40L in T-cells, thereby playing a role in T-cell-dependent antibody responses in activated CD4(+) T-cells and thymus-dependent humoral immunity (PubMed:16936731). Specifically recognizes the MUE3 box, a subset of E-boxes, present in the immunoglobulin enhancer (By similarity). It also binds very well to a USF/MLTF site (By similarity). May regulate lysosomal positioning in response to nutrient deprivation by promoting the expression of PIP4P1 (PubMed:29146937).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).