

Product datasheet for **TL513288**

Lmna Mouse shRNA Plasmid (Locus ID 16905)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Lmna Mouse shRNA Plasmid (Locus ID 16905)
Locus ID:	16905
Synonyms:	Dhe
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Lmna - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 16905). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC015302 , BC094020 , NM_001002011 , NM_001111102 , NM_019390 , NM_001002011.1 , NM_001002011.2 , NM_001002011.3 , NM_019390.1 , NM_019390.2 , NM_019390.3 , NM_001111102.1 , NM_001111102.2
UniProt ID:	P48678
Summary:	This gene encodes a protein that is a member of the lamin family. Nuclear lamins, intermediate filament-like proteins, are the major components of the nuclear lamina, a protein meshwork associated with the inner nuclear membrane. This meshwork is thought to maintain the integrity of the nuclear envelope, participate in chromatin organization, and regulate gene transcription. Vertebrate lamins consist of two types, A and B. This protein is an A-type and is proposed to be developmentally regulated. In mouse deficiency of this gene is associated with muscular dystrophy. Mouse lines with different mutations in this gene serve as pathophysiological models for several human laminopathies. In humans, mutations in this gene lead to several diseases: Emery-Dreifuss muscular dystrophy, familial partial lipodystrophy, limb girdle muscular dystrophy, dilated cardiomyopathy, Charcot-Marie-Tooth disease, and Hutchinson-Gilford progeria syndrome. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, May 2013]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).