

Product datasheet for **TL513396**

Wnt7a Mouse shRNA Plasmid (Locus ID 22421)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Wnt7a Mouse shRNA Plasmid (Locus ID 22421)
Locus ID:	22421
Synonyms:	AI849442; px; tw; Wnt-7a
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Wnt7a - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 22421). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC049093 , BC057586 , NM_009527 , NM_009527.1 , NM_009527.2 , NM_009527.3 , NM_001363757 , NM_009527.4
UniProt ID:	P24383
Summary:	Ligand for members of the frizzled family of seven transmembrane receptors that functions in the canonical Wnt/beta-catenin signaling pathway (PubMed:18230341, PubMed:20530549, PubMed:23629626). Plays an important role in embryonic development, including dorsal versus ventral patterning during limb development, skeleton development and urogenital tract development (PubMed:7885472, PubMed:9769174, PubMed:9790192). Required for central nervous system (CNS) angiogenesis and blood-brain barrier regulation (PubMed:28803732). Required for normal, sexually dimorphic development of the Mullerian ducts, and for normal fertility in both sexes (PubMed:9790192). Required for normal neural stem cell proliferation in the hippocampus dentate gyrus (PubMed:23629626). Required for normal progress through the cell cycle in neural progenitor cells, for self-renewal of neural stem cells, and for normal neuronal differentiation and maturation (PubMed:23629626). Promotes formation of synapses via its interaction with FZD5 (PubMed:20530549). [UniProtKB/Swiss-Prot Function]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).