

## Product datasheet for **TL513440**

### Dcun1d2 Mouse shRNA Plasmid (Locus ID 102323)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Dcun1d2 Mouse shRNA Plasmid (Locus ID 102323)
Locus ID:	102323
Synonyms:	AW495713; mFLJ10704
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Dcun1d2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 102323). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001024504</a> , <a href="#">NM_001042649</a> , <a href="#">NM_001042650</a> , <a href="#">NM_001042651</a> , <a href="#">NM_001042651.1</a> , <a href="#">NM_001024504.1</a> , <a href="#">NM_001024504.2</a> , <a href="#">NM_001042649.1</a> , <a href="#">NM_001042650.1</a> , <a href="#">BC030839</a> , <a href="#">BC052676</a> , <a href="#">BC156825</a>
UniProt ID:	<a href="#">Q8BZJ7</a>
Summary:	Potently stimulates the neddylation of cullin components of SCF-type E3 ubiquitin ligase complexes from the NEDD8-conjugating E2 enzyme UBC12. Neddylation of cullins play an essential role in the regulation of SCF-type complexes activity (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).