

Product datasheet for TL513468V

Bdnf Mouse shRNA Lentiviral Particle (Locus ID 12064)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Bdnf Mouse shRNA Lentiviral Particle (Locus ID 12064)
Locus ID:	12064
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Bdnf - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10^7 TU/ml.
RefSeq:	BC034862, NM 001048139, NM 001048141, NM 001048142, NM 001285416, NM 001285417, NM 001285418, NM 001285419, NM 001285420, NM 001285421, NM 001285422, NM 007540, NM 001048141.1, NM 007540.1, NM 007540.2, NM 007540.3, NM 007540.4, NM 001048139.1, NM 001285416.1, NM 001285417.1, NM 001285418.1, NM 001285419.1, NM 001285420.1, NM 001285421.1, NM 001285422.1
UniProt ID:	<u>P21237</u>
Summary:	The protein encoded by this gene is a member of the nerve growth factor family. It is involved in the growth, differentiation and survival of specific types of developing neurons both in the central nervous system (CNS) and the peripheral nervous system. It is also involved in regulating synaptic plasticity in the CNS. Expression of a similar gene in human is reduced in both Alzheimer's and Huntington disease patients. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar processing to generate mature protein. [provided by RefSeq, Oct 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u> . If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u> .



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GRIGENE Bdnf Mouse shRNA Lentiviral Particle (Locus ID 12064) – TL513468V

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US