

Product datasheet for **TL513704V**

Zhx2 Mouse shRNA Lentiviral Particle (Locus ID 387609)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Zhx2 Mouse shRNA Lentiviral Particle (Locus ID 387609)
Locus ID:	387609
Synonyms:	Afr-1; Afr1; mKIAA0854; Raf
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Zhx2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC059178 , NM_199449 , NM_199449.1 , NM_199449.2
UniProt ID:	Q8C0C0
Summary:	Acts as a transcriptional repressor (PubMed:19515908). Represses the promoter activity of the CDC25C gene stimulated by NFYA (By similarity). May play a role in retinal development where it regulates the composition of bipolar cell populations, by promoting differentiation of bipolar OFF-type cells (PubMed:30146259). In the brain, may promote maintenance and suppress differentiation of neural progenitor cells in the developing cortex (PubMed:19515908).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).