

Product datasheet for **TL513712**

Mc3r Mouse shRNA Plasmid (Locus ID 17201)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Mc3r Mouse shRNA Plasmid (Locus ID 17201)
Locus ID:	17201
Synonyms:	MC3-R
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mc3r - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 17201). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC103669 , BC104730 , BC105667 , BC105668 , NM_008561 , NM_008561.1 , NM_008561.2 , NM_008561.3
UniProt ID:	P33033
Summary:	This gene encodes a member of the melanocortin receptor family. Melanocortin receptors are transmembrane G-protein coupled receptors, which respond to small peptide hormones and exhibit diverse functions and tissue type localization. As part of the central nervous melanocortin system, the encoded protein is competitively bound by either melanocyte stimulating hormone or agouti-related protein to regulate energy homeostasis and adaptation to food restriction. Disruption of this gene results in an increased ratio of weight gain to food intake, increased fat mass, and decreased lean mass, without having a large effect on insulin sensitivity or glucose metabolism. [provided by RefSeq, Dec 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).