

## Product datasheet for **TL513726V**

### Nf2 Mouse shRNA Lentiviral Particle (Locus ID 18016)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Nf2 Mouse shRNA Lentiviral Particle (Locus ID 18016)
Locus ID:	18016
Synonyms:	merlin
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Nf2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC005442</a> , <a href="#">NM_001252250</a> , <a href="#">NM_001252251</a> , <a href="#">NM_001252252</a> , <a href="#">NM_001252253</a> , <a href="#">NM_010898</a> , <a href="#">NM_010898.1</a> , <a href="#">NM_010898.2</a> , <a href="#">NM_010898.3</a> , <a href="#">NM_010898.4</a> , <a href="#">NM_001252253.1</a> , <a href="#">NM_001252252.1</a> , <a href="#">NM_001252250.1</a> , <a href="#">NM_001252251.1</a> , <a href="#">NM_001361675</a> , <a href="#">NM_001361676</a> , <a href="#">NM_001361677</a>
UniProt ID:	<a href="#">P46662</a>
Summary:	Probable regulator of the Hippo/SWH (Sav/Wts/Hpo) signaling pathway, a signaling pathway that plays a pivotal role in tumor suppression by restricting proliferation and promoting apoptosis. Along with WWC1 can synergistically induce the phosphorylation of LATS1 and LATS2 and can probably function in the regulation of the Hippo/SWH (Sav/Wts/Hpo) signaling pathway. May act as a membrane stabilizing protein. May inhibit PI3 kinase by binding to AGAP2 and impairing its stimulating activity. Suppresses cell proliferation and tumorigenesis by inhibiting the CUL4A-RBX1-DDB1-VprBP/DCAF1 E3 ubiquitin-protein ligase complex (By similarity). Plays a role in lens development and is required for complete fiber cell terminal differentiation, maintenance of cell polarity and separation of the lens vesicle from the corneal epithelium.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).