

## Product datasheet for **TL514031**

### Cela2a Mouse shRNA Plasmid (Locus ID 13706)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cela2a Mouse shRNA Plasmid (Locus ID 13706)
Locus ID:	13706
Synonyms:	El; Ela; Ela-2; Ela2; Ela2a
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cela2a - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 13706). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC026552</a> , <a href="#">NM_007919</a> , <a href="#">NM_007919.1</a> , <a href="#">NM_007919.2</a> , <a href="#">NM_007919.3</a>
UniProt ID:	<a href="#">P05208</a>
Summary:	This gene encodes a serine protease enzyme that hydrolyzes elastin. This gene is highly expressed in the pancreatic acinar cells where the encoded preproprotein undergoes processing including signal peptide cleavage to generate an inactive zymogen. The removal of N-terminal activation peptide from the zymogen by trypsin generates active elastase enzyme. This gene is also expressed in the mouse epidermis where it participates in pro-filaggrin processing. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).