

## Product datasheet for **TL514056**

### Arfrp1 Mouse shRNA Plasmid (Locus ID 76688)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Arfrp1 Mouse shRNA Plasmid (Locus ID 76688)
Locus ID:	76688
Synonyms:	1500006I01Rik; AI480700
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Arfrp1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 76688). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC010713</a> , <a href="#">BC046782</a> , <a href="#">BC046786</a> , <a href="#">NM_001165991</a> , <a href="#">NM_001165992</a> , <a href="#">NM_001165995</a> , <a href="#">NM_029702</a> , <a href="#">NR_028585</a> , <a href="#">NM_001165991.1</a> , <a href="#">NM_001165992.1</a> , <a href="#">NM_001165995.1</a> , <a href="#">NM_029702.1</a> , <a href="#">NM_029702.2</a> , <a href="#">NM_029702.3</a> , <a href="#">NM_029702.4</a>
UniProt ID:	<a href="#">Q8BXL7</a>
Summary:	The gene encodes a membrane-associated GTPase that is related to the ADP-ribosylation factor (ARF) and ARF-like (ARL) genes. It plays an essential role in Golgi function controlling recruitment of GRIP domain proteins and ARL1 to the trans-Golgi and trans-Golgi to plasma membrane trafficking of cell surface proteins such as E-cadherin. Deletion of this gene in mice leads to embryonic lethality during early gastrulation, which is at least partly caused by the disruption of E-cadherin trafficking to the cell surface and therefore lack of sufficient cell-cell adhesion in the embryo. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).