

## Product datasheet for **TL514138V**

### **Pcgf2 Mouse shRNA Lentiviral Particle (Locus ID 22658)**

#### **Product data:**

<b>Product Type:</b>	shRNA Lentiviral Particles
<b>Product Name:</b>	Pcgf2 Mouse shRNA Lentiviral Particle (Locus ID 22658)
<b>Locus ID:</b>	22658
<b>Synonyms:</b>	mel-18; Mel18; Rnf110; Zfp144
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>Format:</b>	Lentiviral particles
<b>Components:</b>	Pcgf2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
<b>RefSeq:</b>	<a href="#">BC016419</a> , <a href="#">NM_001163307</a> , <a href="#">NM_001163308</a> , <a href="#">NM_009545</a> , <a href="#">NM_001163307.1</a> , <a href="#">NM_009545.1</a> , <a href="#">NM_009545.2</a> , <a href="#">NM_001163308.1</a>
<b>UniProt ID:</b>	<a href="#">P23798</a>
<b>Summary:</b>	Transcriptional repressor (PubMed:8521824). Binds specifically to the DNA sequence 5'-GACTNGACT-3' (PubMed:8521824). Has tumor suppressor activity (PubMed:8521824). May play a role in control of cell proliferation and/or neural cell development (Probable). Regulates proliferation of early T progenitor cells by maintaining expression of HES1(PubMed:15728456). Also plays a role in antero-posterior specification of the axial skeleton and negative regulation of the self-renewal activity of hematopoietic stem cells (PubMed:8625838, PubMed:15183898). Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility. Within the PRC1-like complex, regulates RNF2 ubiquitin ligase activity (By similarity).[UniProtKB/Swiss-Prot Function]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).