

## Product datasheet for **TL514147**

### Mreg Mouse shRNA Plasmid (Locus ID 381269)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Mreg Mouse shRNA Plasmid (Locus ID 381269)
Locus ID:	381269
Synonyms:	dsu; Gm974; Wdt2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mreg - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 381269). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC068125</a> , <a href="#">NM_001005423</a> , <a href="#">NM_001357488</a> , <a href="#">NM_001005423.1</a> , <a href="#">NM_001005423.2</a>
UniProt ID:	<a href="#">Q6NVG5</a>



[View online »](#)

<b>Summary:</b>	<p>Probably functions as cargo-recognition protein that couples cytoplasmic vesicles to the transport machinery (PubMed:22940130, PubMed:22275436, PubMed:30174147). Plays a role in hair pigmentation, a process that involves shedding of melanosome-containing vesicles from melanocytes, followed by phagocytosis of the melanosome-containing vesicles by keratinocytes (PubMed:15550542, PubMed:3410303, PubMed:22753477). Functions on melanosomes as receptor for RILP and the complex formed by RILP and DCTN1, and thereby contributes to retrograde melanosome transport from the cell periphery to the center (PubMed:22940130, PubMed:22275436). Overexpression causes accumulation of late endosomes and/or lysosomes at the microtubule organising center (MTOC) at the center of the cell (PubMed:19240024, PubMed:30174147). Probably binds cholesterol and requires the presence of cholesterol in membranes to function in microtubule-mediated retrograde organelle transport (PubMed:30174147). Binds phosphatidylinositol 3-phosphate, phosphatidylinositol 4-phosphate, phosphatidylinositol 5-phosphate and phosphatidylinositol 3,5-bisphosphate, but not phosphatidylinositol 3,4-bisphosphate or phosphatidylinositol 4,5-bisphosphate (PubMed:19240024). Required for normal phagosome clearing and normal activation of lysosomal enzymes in lysosomes from retinal pigment epithelium cells (PubMed:19240024). Required for normal degradation of the lipofuscin component N-retinylidene-N-retinylethanolamine (A2E) in the eye (PubMed:19240024). May function in membrane fusion and regulate the biogenesis of disk membranes of photoreceptor rod cells (Probable).[UniProtKB/Swiss-Prot Function]</p>
<b>shRNA Design:</b>	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a>.</p>
<b>Performance Guaranteed:</b>	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>