

## **Product datasheet for TL514320**

### **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

# Ctnna2 Mouse shRNA Plasmid (Locus ID 12386)

**Product data:** 

**Product Type:** shRNA Plasmids

Product Name: Ctnna2 Mouse shRNA Plasmid (Locus ID 12386)

**Locus ID:** 12386

Synonyms: Al481747; Catna; Catna2; cdf; chp

**Vector:** pGFP-C-shLenti (TR30023)

E. coli Selection: Chloramphenicol (34 ug/ml)

**Mammalian Cell** 

Selection:

Puromycin

Format: Lentiviral plasmids

Components: Ctnna2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 12386).

5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.

**RefSeq:** BC079648, NM 001109764, NM 009819, NM 145732, NM 001355193, NM 001109764.1,

NM 009819.1, NM 009819.2, NM 145732.2, BC051995

UniProt ID: 061301

**Summary:** May function as a linker between cadherin adhesion receptors and the cytoskeleton to

regulate cell-cell adhesion and differentiation in the nervous system (PubMed:12123610, PubMed:15034585). Required for proper regulation of cortical neuronal migration and neurite growth. It acts as negative regulator of Arp2/3 complex activity and Arp2/3-mediated actin polymerization. It thereby suppresses excessive actin branching which would impair neurite growth and stability (By similarity). Regulates morphological plasticity of synapses and cerebellar and hippocampal lamination during development. Functions in the control of

startle modulation (PubMed:12089526).[UniProtKB/Swiss-Prot Function]

**shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u>.





# Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).